



ENTONOX

MIDOGAS

We make a range of Patient Breathing Circuits to suit our MIDOGAS Analgesic Units as well as a circuit suitable for use with the Entonox System.

Our range includes:

- Standard Breathing Circuit (single BLUE breathing hose)
- Breathing Circuit with Scavenge Hose (double circuit with BLUE and PINK hoses)
- Entonox Circuit (single PINK scavenge hose)

Our breathing circuits feature:

- Two fitted one-way valves to eliminate the risk of breathing exhaled air
- Ergonomically designed hand-piece for easy grip, stability and comfort
- Streamlined parallel circuit hoses

As we manufacture all our breathing circuits, we would be happy to discuss any specific requirements you may have and customise a circuit to meet your needs.

**Easy Grip and
Comfortable Hold**



Disposable Breathing Circuits

Maternity Delivery Suite & Emergency

Great Aussie Made Products
with local service & technical support



www.pdmedical.com.au



PDMedical Pty Ltd
80-84 Arkwright Drive,
Dandenong South,
Victoria, Australia.
T: +61 3 9708 9708

E: sales@pdmedical.com.au

Part of the
MIDOGAS
Range of Products
Made in Australia 

Expandable tubing circuits can be configured to your specific requirements

