

## Minimise Pathogen Transmission

normal - airborne - contact - droplet

The Personal Protective Equipment (PPE) Caddy has been developed in consultation with IC nurses in hospitals around Australia to help nurses, staff and visitors comply with safety precautions in standard and isolation scenarios in order to improve safety and minimise the risk of pathogen transmission.

It holds all the PPE necessary for **normal, airborne, contact** and **droplet** precautions plus **donning signage** for immediate patient care in a normal or isolation environments.

Complies with the Australian Safety and Quality Standard on PPE, signage and CDC infection control protocols.

Declutters hallways of hospital and other healthcare facilities and provides a consistent, clean and professional layout.

Available in white or soft grey as a wall unit or as a multi-functional mobile unit.

Easy to clean and keep clean

Made by us in Australia and distributed direct so it is always available.



# PPE Caddy

HOSPITALS  
AGED CARE FACILITIES  
GP CLINICS

Great Aussie Made Products  
with local service & technical support



Innovative : Helpful : Trusted  
[www.pdmedical.com.au](http://www.pdmedical.com.au)



PDMedical Pty Ltd  
80-84 Arkwright Drive,  
Dandenong South, 3175  
Victoria, Australia.  
T: +61 3 9708 9708  
E: [peter@pdmedical.com.au](mailto:peter@pdmedical.com.au)

Accommodates multiple suppliers' PPE products so it eliminates the need for maintenance to remove/ replace specific holders and fixings and patch up damaged walls

## Glide your PPE to the Point of Care with our Multi-Functional Mobile PPE Unit



The Mobile PPE Unit is compact and designed to hold all the PPE on the **donning side** necessary for normal, airborne, contact and droplet precautions as well as a clipboard to hold the appropriate signage.

The **procedure side** of the mobile unit contains compartments to cater for:

- vomit bags,
- wipes,
- sharps containers and binders.

It also contains a removable tray for use by nurses and convenience when carrying out a procedure.

It also holds the doffing sign on the procedure side as a constant reminder for compliance.

The mobile unit has handle either side to decrease cross contamination.

The base has a low centre of gravity for safety and stability and six easy roll rubber castors.

## Multi-Functional Mobile PPE Unit

## Ideal for Hospital Wards and Aged Care Facilities

**Great Aussie Made Products**  
with local service & technical support



**Innovative : Helpful : Trusted**  
[www.pdmedical.com.au](http://www.pdmedical.com.au)



PDMedical Pty Ltd  
80-84 Arkwright Drive,  
Dandenong South, 3175  
Victoria, Australia.  
T: +61 3 9708 9708  
E: [peter@pdmedical.com.au](mailto:peter@pdmedical.com.au)

The ideal mobile unit to help the nurses do their job, comply with your protocols to provide better and safer care.

### Simple, Intuitive, Instructive

Easy step by step instructions help visitors and staff comply



Where a full gown is required a wire basket can be clipped on to the bottom of the caddy

We also have a **Clean-Up Caddy**



It also fits tubs with side access to wipes



The Personal Protective Equipment (PPE) Caddy is designed to hold all the PPE required for normal and other levels of isolation.

The arrangement of the PPE in the caddy is in order of the donning sequence prescribed by the CDC and Australian Safety Standards.

The PPE items have been carefully placed in compartments within the caddy so that the donning sequence follows a clockwise direction from bottom left to right.

These compartments are numbered 1 to 5 with a soft emboss and correlate to the numbers on the donning signage


This helps in training new and agency nurses and serves as memory aid for all nurses to comply with protocols.

Easy to demonstrate, easy to understand, easy to remember, easy to comply.

## PPE Compartments

If you'd like some training material to help you help your nurses and staff please contact us. We'd be happy to help.

Our unique signs are available ready for you to print or we can customise and supply them ready printed on hard wearing plastic.



Adhering to the right protocols in the various isolation cases is pivotal to the minimisation of pathogen transmission within a room, a ward and medical facility.

Many hospitals have agency nurses, casual staff, general healthcare workers, and doctors who work in other hospitals which may have different protocols and expectations.

Ensuring that your facility's protocols and expectations are met by visitors and staff alike requires clear signage and easily accessible PPE.

The PPE Caddy incorporates a clipboard arrangement for the donning and doffing instruction sheet suited to the situation to be posted at the point of use.

Signs can be customised to reflect your requirements and protocols.



The donning sequence numbering (1-5) on the signs corresponds with the embossed numbers on the PPE caddy (1-5) to maximise clarity and ensure full compliance.



## PPE Caddy Donning/Doffing Signage

Signs are easily interchangeable so you can react quickly and uniformly.

### Components



PPE Caddy



Poster Clipboard



"Hook-On" Basket  
(for Full Gowns)



"Clean-Up"  
Caddy

The PPE Caddy can be configured in a number of ways to suite your requirements.

#### PPE Caddy - Wall mounted

- PPE Caddy
- Interchangeable poster clipboard
- Basket for gowns

#### Clean-Up Caddy

- Vomit bags
- Wipes  
*(holds most commonly available wipes in satchels and tubs)*

### Assembly Options



"Easily hooked on  
and off"



#### Mobile PPE Unit

- PPE Caddy
- Interchangeable poster clipboard
- Mobile stand
- Gown Basket
- Cleanup Caddy:
  1. Vomit bags
  2. Wipes compartment
- Bracket to fit:
  1. Sharps container
  2. Folder
- Dispensing tray

All units accommodate a wide variety of widely used PPE equipment, from hand wash bottles, wipes, apron, masks, eyewear and gloves.

## PPE Caddy Options

We can customise an assembly to satisfy your specific needs.  
Contact us now.




Fixed wall mounted caddy  
installs with four screw  
fixtures.

**PPE Caddy:**  
Size (cm):  
60H x 39W x 10D



Mobile unit is supplied  
sub-assembled

**Mobile Stand:**  
Size (cm):  
60D x 45W x 150H



**"Clean-Up" Caddy**  
Size (cm):  
15H x 39W x 10D



**"Hook On" Basket**  
Size (cm):  
20H x 39W x 10D

The Personal Protective Equipment (PPE) Caddy is made from hard wearing, lightweight and strong plastic.

The PPE unit is easily installed on a wall or door as is the Clean Up Caddy.

The both the PPE and Clean-Up caddys accommodate multiple suppliers' PPE product so it eliminates the need for maintenance to remove/replace specific holders and fixings and patch up damaged walls.

The full gown basket can easily be hooked on to the bottom of the PPE or Clean up caddy when needed. The baskets stack well and only need a small area for storage.

The mobile unit is supplied sub-assembled and requires no special skillsets to assemble with simple hand tool.

This universal PPE Caddy accommodates multiple PPE supplier products which minimises damage to your walls from constant changes to holders.

Keeps maintenance happy too.

PPE Caddy Installation

Fix and forget