

Designed to  
clear the air &  
protect your staff.

Our Gas Scavenge Unit's scavenge buffer chamber enables the hospital's "in wall" suction or scavenge facilities to safely draw the expired analgesic gases from a patient's breathing circuit in a safe and simple manner, thereby reducing work place exposure to Nitrous Oxide.

Provides a safe and cost effective solution for hospital environments, creating a safer workplace for your staff.

The Scavenge Unit can be wall-mounted next to MIDOGAS Units if required or incorporated as part of our Mobile Stand.



## MIDOGAS Gas Scavenge Unit

Simple  
Effective  
Integrated

**Great Aussie Made Products**  
with local service & technical support



[www.pdmedical.com.au](http://www.pdmedical.com.au)



PDMedical Pty Ltd  
80-84 Arkwright Drive,  
Dandenong South,  
Victoria, Australia.  
T: +61 3 9708 9708  
F: +61 3 9793 8084  
E: [peter@pdmedical.com.au](mailto:peter@pdmedical.com.au)

Part of the  
**MIDOGAS**  
Range of Products  
Made in Australia 