

Protect your nurses
from contaminated
sharps and
needlestick injuries



The Sharps Caddy is the safest way to load, carry and administer drugs and dispose of contaminated sharps and needles at the point of use.

Ergonomically friendly and easy to clean, it provides a safe and organised workstation.


The Sharps Caddy includes an injection tray and compartments to hold a sharps container and a variety of hand hygiene bottles.

- Safer for clinicians and patients
- Easy to clean
- Strong, compact, stable
- Perfectly organised mobile workstation

The A.T. Sharps Caddy with the larger tray is

ideally suited for aseptic procedures including:

- Peripheral IV cannulation
- Blood collection
- Administering IV medication
- Injecting intra-muscular / sub-cutaneous



SHARPS CADDY

Safer
Cleaner
Organised

Great Aussie Made Products
with local service & technical support



www.pdmedical.com.au



PD Medical Pty Ltd
80-84 Arkwright Drive,
Dandenong South,
Victoria, Australia.
T: +61 3 9708 9708
E: sales@pdmedical.com.au

Promotes hand
hygiene
compliance

Reducing Sharps Injuries

In the Hospital Ward



... The Sharps Caddy

A caddy for use in administering drugs and the safe carry and **disposal** of contaminated sharps **at the point of use**.



Disposal of sharps at the point of use makes the Sharps Caddy effective in protecting nurses from the risk of injury and infection.

A simple solution to a difficult problem



Nurse love them...

"It's intuitive to use and having the sharps bin and hand wash next to the tray allows us to do what comes naturally. A constant reminder of best safe practice with everything you need right next to you!"

SHARPS CADDY

Safer
Cleaner
Organised

Great Aussie Made Products
with local service & technical support



www.pdmedical.com.au



PDMedical Pty Ltd
80-84 Arkwright Drive,
Dandenong South,
Victoria, Australia.
T: +61 3 9708 9708
E: sales@pdmedical.com.au


Promotes hand
hygiene
compliance

A simple solution to 2 difficult problems

Needle stick injury prevention & Hand hygiene compliance

We understand the pressure nurses are under, so we engineered the “Sharps Caddy” to simplify the patient care process and provide convenient access to hand rub and sharps container to help nurses:-

- * do their job quicker,
- * dispose of the contaminated needle immediately following use.
- * clean hands immediately before the medication is administered as well as immediately after



The “Sharps Caddy” provides a complete ergonomic solution **at the point of care**.

Designed, developed and made in Australia



The **Compact Sharps Caddy** can be used for a variety of procedures including:

- Intravenous antibiotics
- Subcutaneous Injections
- Intramuscular Injections



The **A.T. Sharps Caddy** with the larger tray is ideally suited for aseptic procedures including:

- Peripheral IV cannulation
- Blood collection
- Administering IV medication
- Injecting intra-muscular / sub-cutaneous



Contact: peter@pdmedical.com.au



Sharps Caddy



Caddy: SC-100B-STV

Caddy: SC-100B-STV

Caddy: SC-100P-STV

Sharps Container:
SC-CONT-1.4LT

Cytotoxic Container:
CC-CONT-1.4LT

Sharps Container:
SC-CONT-1.4LT

Insert: SC-INS-200B

Caddy: SC-AT-200B

Item (P=Pink, B=Blue, PP=Purple)	Qty/Box	
Small Sharps Caddy		
Code: SC-100B-STV	4	
Code: SC-100P-STV	4	
Code: SC-100PP-STV	4	
AT Sharps Caddy		
Code: SC-AT-200B	4	
Code: SC-AT-200P	4	
Code: SC-AT-200PP	4	
Sharps Containers		
Code: SC-CONT-1.4LT	30	Sharps
Code: CC-CONT-1.4LT	30	Cytotoxic
Inserts		
Code: SC-INS-200B	To match # Caddys	
Trial Pack		
2 x Caddys (Small or Large)		
+ 6 x Containers		

Our Sharps Caddy is designed to fit most of the common brands of hand-wash or hand sanitizer.

These bottles fit nicely into the caddy pocket:

- 3M Avagard Antiseptic Handrub - 500mL
- Cutan Alcohol Foam - Bottle 400mL
- Purell Antiseptic Hand Gel - 350mL

You will need a caddy insert if you use the following hand wash bottles:

- Aqium Gel Antibacterial Hand Gel - 375mL
- Microshield Angel Blue - 500mL

Note: Bactol Hand Sanitizer 500mL does not fit the caddy.

Please be sure to advise us of the hand-wash bottle you use so we can determine if an insert is required

