

Smooth bored
to reduce
turbulence




Our tapered Adaptors and Connectors (straight and Y) are manufactured from clear polypropylene with clear stepped ring terminations to provide strong, secure connections.

They are smooth bored with a gentle taper to minimise turbulence.

Supplied as individually packed sterile products or non-sterile packs.

Our Tapered Spigots are manufactured from polypropylene. They provide a versatile method of plugging a wide variety of sizes of tubing terminations.

Manufactured in Melbourne so we'd be happy to discuss and tailor a pack specific to your needs.



Tube Connectors, Adaptors & Spigots

High quality & wide range of sizes

Great Aussie Made Products
with local service & technical support



www.pdmedical.com.au







PDMedical Pty Ltd
80-84 Arkwright Drive,
Dandenong South,
Victoria, Australia.
T: +61 3 9708 9708

E: peter@pdmedical.com.au

Single pack sterilised or loose packed in bulk

Tube Connectors, Adaptors & Spigots

High quality &
wide range
of sizes

TUBE CONNECTOR: Small : 4-7mm OD x 3mmID	
TUBE CONNECTOR: Medium : 7-10mm OD x 6mmID	
TUBE CONNECTOR: Large : 10-14mm OD x 8mmID	
TUBE ADAPTOR: X-Small/Medium : 3-5mm/7-10mmOD X 2mmID	
TUBE ADAPTOR: Small/Medium : 4-7mm/7-10mmOD X 3mmID	
TUBE ADAPTOR: Small/Large : 4-7mm/10-14mmOD X 3mmID	
TUBE ADAPTOR: Medium/Large : 7-10mm/10-14mmOD X 6mmID	
Y-TUBE CONNECTOR: Small : 5-8mmOD X 4mmID	
Y-TUBE CONNECTOR: Medium : 8-10mmOD X 6mmID	
Y-TUBE CONNECTOR: Large : 12-14mmOD X 10mmID	
SPIGOT 10mm : 0-10mmOD x 49mmLong	
SPIGOT 13mm : 0-13mmOD x 52mmLong	